

# "I've stopped falling!"

## Finally! a walking aid for Parkinson's

### Designed for Parkinson's Disease!

#### • Stability

We created the *U-Step Walking Stabilizer* to increase your independence. The secret is in our patented U-shaped base. Its ultra-stable foundation braces you in every direction. It is not like pushing a walker. Instead, the *U-Step* surrounds you and moves with you. You will feel as stable as you would feel holding onto another person's arm.

#### • Safety

The innovative braking system is easy to use and puts you in complete control. The *U-Step* will not roll unless you are ready to walk. When you lightly squeeze a hand brake, the unit will roll with you. Once you release the hand brake the unit will stop immediately. This feature is particularly helpful when standing up from a chair because, the unit will not roll away from you.

#### • Rolls easily

#### • Improves posture

#### • Improves walking stride

#### • Rolling resistance control

#### • Folds for transport

#### • Easy turning



Optional feature:  
red button for operating  
LaserLight™

### Laser Light™ helps prevent freezing

If you suffer from freezing episodes, *U-Step* has a solution. Introducing *LaserLight™* an exclusive optional feature of the *U-Step*. Simply press the red button on the handlebar and a bright red laser line is projected on the floor for you to step over. See projected line illustrated on back page. It is amazing! Our laser offers an entirely safe, obstacle-free, visual cue that helps you break the "freezing episode," and walk normally with increased stride length.

*"My father shows off the laser to all his cohorts with Parkinson's. He says that if you can throw down a pencil and step over it, the laser will work for you. He does not use the urinal at the side of his bed in the middle of the night any more. With the U-Step and the laser he can make it to the washroom."*

— Dr. Carlo Klott, M.D., Austin, TX

Note: None of the physicians or therapists quoted in this material have any financial interest in selling the U-Step.

**Stop falling! Order your U-Step and LaserCane NOW!**

**LaserCane is ideal for those seeking minor walking support**

# All the features you wanted in an advanced walker!



## • Tension Adjustment

By turning the knobs at the base of the unit, you can regulate your walking speed. This feature is particularly comforting if you are afraid of a rolling walker "getting away" from you.

## U-Step Walking Stabiliser dimensions

Weight capacity	170	Kg
Height adjustment range	155-185	cm
Overall width	58	cm
Width of large padded seat	51	cm
Height from floor to seat	56	cm
Length	64	cm
Turning circle	74	cm
Weight	8.6	kg
Folded dimensions (approx.)	91x57x 41	cm
Material of frame	tubular steel	
Colour choices	black or blue	



*"The braking system is amazing and the laser line gets his coordination back. [The laser] gets him to lift his feet; he drags his feet without the light."*

— Mrs. Robert Erickson, Traverse City, MI

*"Great! Cliff uses the U-Step during his "off" times. He is able to care for himself instead of depending on me for support."*

— Mrs. Norma Eash



Folds for transport.

Made in America!

**△ Recommended at dozens of major Parkinson's centers! Used by thousands of people with Parkinson's**

NSW Distributer of U-Step Walker and Lasercane

**Roland Priestley Pty Ltd**

**26/12 Cecil Rd**

**Hornsby NSW 2077**

**PH: 02 9482 4893**

[www.rolandpriestley.com.au](http://www.rolandpriestley.com.au)